

ODIN WORKSHOP · STRASBOURG · 27 SEPTEMBER 2026

Current and Future Applications of **XR** and **AI** in Dentistry



Roberto Mangano

MSc Computer Engineering
Digital Dentistry Society Certified Speaker

www.digital-dentistry.org



Digital
Dentistry
Society

Virtual Reality

Augmented Reality

AI Glasses

Mixed Reality

TODAY'S QUESTIONS

What can we do today with these technologies?

How to navigate the current market?

Where does AI fit into it?



Current and Future Applications of XR and AI in Dentistry

01

Key concepts

- AR / VR / MR
- How to navigate the current market

02

Current applications

- VR Training
- VR Simulation
- Virtual Desktop
- University education

03

Mixed Reality in Dentistry

- The Hololmplant case study

04

Dynamic Navigation with Mixed Reality

- In vitro and in vivo studies




REALITY-VIRTUALITY CONTINUUM


REALITY



AUGMENTED REALITY (AR)



MIXED REALITY (MR)



VIRTUAL REALITY (VR)



REALITY

EXTENDED REALITY CONTINUUM

VIRTUALITY




REALITY

What we see and experience naturally



AUGMENTED REALITY

Adds digital information, images, text on top of the real world



MIXED REALITY

Digital objects exist and interact with the real world in real time




VIRTUAL REALITY

A complete immersive virtual environment with no real-world visuals



AR vs. VR: Key Differences



AR Augmented Reality

- Adds digital elements to a live view through a smartphone or AR glasses
- Interactive: overlays information onto the real world



VR Virtual Reality

- Immerses the user in a fully virtual environment, cut off from the real world
- Used for simulation, training and gaming



AR: See-through glasses

This technology allows users to see the real world directly through a transparent lens, onto which digital content is overlaid.

EXAMPLES

- HoloLens 2, Magic Leap 2, Xreal

ADVANTAGES

- Natural view of reality
- Lightweight and comfortable
- Limited brightness of virtual elements



Xreal



HoloLens 2



Magic Leap 2



AR/AI Glasses

Integrate AI for real-time understanding of the environment (pictures, videos), contextual assistance, and natural interaction through voice, vision and sensors.

EXAMPLES

- Ray-Ban Meta Smart Glasses, Alibaba Qwen, HTC Vive Eagle, Meta Display

FEATURES

- Computer Vision & AI integration (Gemini, OpenAI, Meta AI)
- Audio answers, automatic translation, music playback
- Always require internet connectivity
- Newer generations feature a display



Meta Ray-Ban



Meta Display AI Glasses



HTC Vive Eagle



TCL X2



AR: Pass-through headsets

Pass-through technology uses **external cameras** to capture the real world, which is then displayed on **opaque screens** inside a VR headset.

EXAMPLES


- Meta Quest, HTC Vive, Apple Vision Pro

ADVANTAGES

- Full control over displayed content
- Better integration of virtual and real
- Versatility and better contrast



Meta Quest 3



Apple Vision Pro



Meta VR



See-through vs. pass-through

APPLICATION SCENARIOS



See-through

Preferred where the real world must be seen clearly, without distortion

BEST FOR

- **Medical AR**
- **Military**
- **Industrial applications**




Pass-through

Ideal for mixed reality that needs a seamless transition between AR and VR

BEST FOR

- **Training and simulations**
- **Immersive gaming**
- **Creative design**





Current and Future Applications of XR and AI in Dentistry

01

Key concepts

- AR / VR / MR
- How to navigate the current market

02

Current applications

- **VR Training**
- **VR Simulation**
- **Virtual Desktop**
- **University education**

03

Mixed Reality in Dentistry

- The Hololmplant case study


04

Dynamic Navigation with Mixed Reality

- In vitro and in vivo studies








MODULAR VR SYSTEM
VR CLASSIC





MODULAR VR SYSTEM
R&D PROJECT
VR + HAND TRACKING



VR SIMULATION

Simulation with haptic feedback





MIXED REALITY

Intraoral scanner + Virtual Desktop



TINALP

University education ([Link](#))



**Universität
Zürich** UZH

Plattenstrasse



The Digital Dentistry Ecosystem

#SCAN!

Data Acquisition

From real to virtual

- Intraoral Scanners (IOS)
- Face Scanners (FS)
- Cone Beam Computed Tomography (CBCT)
- Mandibular Motion Tracking



#PLAN!

Data Processing

Software & apps

- **Artificial Intelligence (AI)**
- **Virtual / Augmented / Mixed Reality (VR/AR/MR)**



#MILL! #MAKE!

Production

From virtual to real

- Milling units
- 3D Printers
- Selective Laser Melting (SLM)



#DONE!

Clinical Application

Chairside

- **Mixed Reality-based Dynamic Navigation**
- Robotics



Current and Future Applications of XR and AI in Dentistry

01

Key concepts

- AR / VR / MR
- How to navigate the current market

02

Current applications

- VR Training
- VR Simulation
- Virtual Desktop
- University education

03

Mixed Reality in Dentistry

- **The Hololmplant case study**

04

Dynamic Navigation with Mixed Reality

- In vitro and in vivo studies



Our Protocol

- 1 Intraoral Scan (**IOS**) + Cone Beam Computed Tomography (**CBCT**)
- 2 Artificial Intelligence (**AI**) for automated file alignment and segmentation
- 3 Virtual Diagnostic Waxing
- 4 Implant Planning with Mixed Reality (MR)**
- 5 Case Validation Implant Planning Software (2D/3D)
- 6 Design and Fabrication of the Surgical Guide
- 7 Surgery on Patient
- 8 Prosthetic Phase: Intraoral Scan (**IOS**)
- 9 Computer-Assisted-Design (CAD) of Prosthetic Restorations
- 10 Delivery of the Prosthetic Restorations


F Mangano, K R Yang, H Lerner, O Admakin, C Mangano

Artificial Intelligence and Mixed Reality for Dental Implant Planning: a technical note


Clinical Implant Dentistry & Related Research, 2024;26:942–953




Cone Beam Computed Tomography (CBCT)



Relu segmentation & alignment



The standard process



Implant Planning with Mixed Reality (MR)




STEP 2B

The complete set of **STL files** is uploaded in the HoloImplant software



Achieved results and research products in the period



FEATURES

- Multiple rendering modes for DICOM files
- Manage DICOM and 3D files at the same time



Finanziato
dall'Unione europea
NextGenerationEU



Ministero
dell'Università
e della Ricerca




Italiadomani
PIANO NAZIONALE
DI RIPRESA E RESILIENZA



SAPIENZA
UNIVERSITÀ DI ROMA



Development of XR-Optimized Rendering Techniques



FEATURES

- Slicing planes to explore inside DICOM files
- Generate new slices to explore DICOM data
- Optimized AMV3D file for volume rendering



Finanziato
dall'Unione europea
NextGenerationEU



Ministero
dell'Università
e della Ricerca



Italiadomani
PIANO NAZIONALE
DI RIPRESA E RESILIENZA



SAPIENZA
UNIVERSITÀ DI ROMA



Scientific basis for our alignment and segmentation

Our IOS–CBCT pipeline is built on these peer-reviewed methods:

ALIGNMENT

Besl & McKay • IEEE TPAMI 1992

A method for registration of 3-D shapes

Original ICP: our automatic fit step

Chetverikov et al. • Image Vis Comput 2005

Robust Euclidean alignment of 3D point sets: the trimmed ICP algorithm

Trimmed ICP: robust when the scan covers only part of the CT

Horn • JOSA A 1987

Closed-form solution of absolute orientation using unit quaternions

Rigid fit from landmark pairs

Flügge et al. • Clin Oral Implants Res 2017

Registration of CBCT data and intraoral surface scans: a prerequisite for guided implant surgery

Why registration accuracy matters clinically

SEGMENTATION

Pauwels et al. • Dentomaxillofac Radiol 2015

CBCT-based bone quality assessment: are Hounsfield units applicable?




CBCT grey values are not true HU

Pauwels et al. (SEDENTEXCT) • Br J Radiol 2013

Variability of dental cone beam CT grey values for density estimations

Grey values vary across 13 CBCT devices, so one fixed threshold does not transfer





Load scan(s)...

STL, PLY, OBJ(+MTL+textures), GLB — select all files for one scan together.

Case_3_Maxilla_scan.stl (229,885 tris)

Load DICOM file...

Load DICOM folder...

251 slices, 512x512, spacing 0.352/0.352/0.500mm

LAYERS

CT surface (3D view)

Point size

HU threshold
1500 HU — 103,154 points

Window center

Window width

Cut planes (3D)

PICK LANDMARKS (ACTIVE SCAN)

Click a point on the **active scan** (3D view), then the matching point on the **CT surface**. Repeat for at least 3 pairs for a stable fit.

Undo last click
Clear all pairs

PAIRS (0)

No pairs yet.

LANDMARK FIT (REFERENCE)

Pick at least 3 pairs to compute a fit.

MANUAL TRANSFORM (ACTIVE SCAN — THE REAL ANSWER)

Mean distance, scan **1.537 mm (2360/2508 sample pts matched)** — CT surface



State of the art: the Relu process

01

Segment the CBCT

A 3D U-Net segments the jaw bones and labels each tooth individually.

39 s vs 133 min manually^{1,2}

Dice 92.6% (bone)



02

Isolate the crowns

Crowns are segmented before matching, so metal artefacts from restorations are ignored.

0.03 mm surface deviation³

Dice 0.91, heavy artefacts



03

Register IOS onto CBCT

The intraoral scan is fitted automatically onto the segmented teeth.

51 s, median 0.04 mm^{4,5}

7.3× faster than experts

Result: ≈0.03–0.2 mm vs expert registration, in under 1 minute

SOURCES

1. Preda F et al. *Automated segmentation of the maxillofacial complex from CBCT.* **J Dent** 2022

2. Shaheen E et al. *Multi-class tooth segmentation and classification on CBCT.* **J Dent** 2021



3. Elgarba BM et al. *AI crown segmentation for IOS-to-CBCT registration in high-artifact scenarios.* **J Prosthet Dent** 2025

4. Elgarba BM et al. *AI-based automated multimodal registration of CBCT and intraoral scan.* **Clin Oral Implants Res** 2024

5. Preda F et al. *Validation of automated registration of intraoral scan onto CBCT.* **J Dent** 2024



Implant Planning with Mixed Reality (MR)



HOLOIMPLANT

Mixed reality at the chair



What's next?

COVERED TODAY

- 01 Key concepts
AR/VR/MR · How to navigate the current market
- 02 Current applications
- 03 Mixed Reality in Dentistry: The Hololmplant case study

04

NEXT

Dynamic Navigation with Mixed Reality



MR in Implant-Guided Surgery

We have **three different** guided surgery options:

01



Static guided surgery

02



Dynamic navigation

03




**Mixed Reality-based
Dynamic Navigation**



HOLOLENS I TESTING – CUSTOM SOFTWARE

DNS integration (2016)



CASE REPORT

Open Access

Augmented reality for dental implantology: a pilot clinical report of two cases



Gerardo Pellegrino^{1*}, Carlo Mangano², Roberto Mangano³, Agnese Ferri¹, Valerio Taraschi⁴ and Claudio Marchetti⁵



Our Protocol with MR-based Dynamic Navigation

SURGICAL PHASE

- | | |
|--------|---|
| #SCAN! | 1. Intraoral Scan (IOS) + Cone Beam Computed Tomography (CBCT) |
| #PLAN! | 2. Artificial Intelligence (AI) for File Alignment/Segmentation and Implant Planning |
| #DONE! | 3. Surgery on Patient with Mixed Reality-based Dynamic Navigation |

PROSTHETIC PHASE

- | | |
|--------|--|
| #SCAN! | 4. Prosthetic Phase: Intraoral Scan (IOS) |
| #PLAN! | 5. Computer-Assisted-Design (CAD) of Prosthetic Restorations |
| #MILL! | 6. Manufacturing of the Restorations |
| #DONE! | 7. Delivery of the Prosthetic Restorations |

Elgarba M, Fontanele R, Mangano F, Jacobs R.
Novel AI-based automated virtual implant placement: Artificial versus Human Intelligence
Journal of Dentistry 2024

Shusterman A, Nashef R, Tecco S, Mangano C, Lerner H, Mangano F.
Accuracy of Implant Placement using a Mixed Reality-Based Dynamic Navigation System versus Static Computer-Assisted and Freehand Surgery: an in Vitro Study.
Journal of Dentistry 2024

Shusterman A, Nashef R, Tecco S, Mangano C, Mangano F.
Implant placement using Mixed Reality-Based Dynamic Navigation: a proof of concept.
Journal of Dentistry 2024



MR-DN: *in vitro* study design

STUDY DESIGN

In Vitro Methodology & Study Design

45 Partially Edentulous Maxillary Models

Teeth missing at #15, #16, #25. Preoperative CBCT and Intraoral Scans merged in Implant Planner software.

Group 1 (MR-DN)
15 models




Group 2 (s-CAIS)
15 models

Group 3 (Freehand)
15 models



135 Total Dental Implants Placed

All osteotomies performed by the same experienced clinician to eliminate variable operator skill. 3.75 × 10 mm MIS Lance+ implants used.

THREE SURGICAL TECHNIQUES



MR-DN: *in vitro* setup and results



The Diagnostic Scorecard (Statistical Significance)

	MR-DN vs. Freehand	MR-DN vs. s-CAIS
XY Shift	👑 MR-DN Superior ($P < 0.05$)	🛡️ s-CAIS Superior
Z Depth	👑 MR-DN Superior ($P < 0.05$)	👑 MR-DN Superior
3D Entry	👑 MR-DN Superior ($P < 0.05$)	⚖️ Statistically Equivalent
3D Apex	👑 MR-DN Superior ($P < 0.05$)	👑 MR-DN Superior
Angulation	👑 MR-DN Superior ($P < 0.05$)	👑 MR-DN Superior



The Metrics of Deviation

3D Entry Error (En)


The linear coronal deviation at the implant platform.

3D Apex Error (Ap)

The linear deviation at the most apical point of the implant.

Angulation Error (An)

The 3D angular discrepancy between the planned longitudinal axis and the actual placed axis.




PUBLISHED STUDY

Accuracy of Implant Placement using a Mixed Reality–Based Dynamic Navigation System versus Static Computer–Assisted and Free–hand Surgery: an in Vitro Study.

Shusterman A, Nashef R, Tecco S, Mangano C, Lerner H, Mangano F. [Journal of Dentistry](#) 2024



Average Implant Deviation: The Macroscopic View



STUDY RESULTS

Across all three spatial metrics, **MR-DN dramatically outperformed Freehand placement**, while matching or slightly exceeding the precision of static guides.



The accuracy scorecard

MR-DN vs. Freehand

Is Mixed Reality statistically superior to human estimation?

3D Entry Deviation (En)	Mean difference -0.75mm	 P < 0.0001
3D Apex Deviation (Ap)	Mean difference -1.12mm	 P < 0.0001
Angulation Deviation (An)	Mean difference -3.58°	 P < 0.0001
Z-Axis (Depth)	Mean difference -0.43mm	 P < 0.0001

MR-DN demonstrates **absolute statistical superiority** over the **Freehand method** across every measured parameter, effectively eliminating the human variance associated with analog site preparation.

MR-DN vs. Static Guides


Does Mixed Reality match the rigid control of physical templates?

3D Entry Deviation (En)		 Equivalent (P = 0.6536) No statistical difference.
3D Apex Deviation (Ap)	Mean difference -0.09mm	 MR-DN Superior (P = 0.0896 T-Test / 0.0327 Wilcoxon)
Angulation Deviation (An)	Mean difference -0.35°	 MR-DN Superior (P = 0.0134)
Z-Axis (Depth)	Mean difference -0.20mm	 MR-DN Superior (P < 0.0001)
XY-Axis (Lateral)	Mean difference +0.14mm	 s-CAIS Superior (P = 0.0040)

MR-DN achieves clinical parity with **static guides** at the **platform**, while offering **statistically significant improvements in depth control, apical positioning, and trajectory angle**.



MR-DN: *in vivo* study







PUBLISHED STUDY


Implant placement using Mixed Reality–Based Dynamic Navigation: a proof of concept

Shusterman A, Nashef R, Tecco S, Mangano C, Mangano F. [Journal of Dentistry](#) 2024





MR-DN: a real alternative to guided surgery



MERGED WITH THE REAL ENVIRONMENT

Implant planning

3D models

Computer screens

Calibrated surgical instruments

These are superimposed on the patient, so the surgeon can follow the real-time information from the MR-DN system and the surgical site simultaneously.

The use of Mixed Reality Dynamic Navigation can change the perspectives of guided dental implant surgery as an alternative to the classic static and dynamic guided surgical techniques

Mirela Feraru, Ariel Shusterman, Carlo Mangano, Francesco Guido Mangano

Clinical Integration of Artificial Intelligence and Mixed Reality-Based Dynamic Navigation for Implant Placement in the Anterior Maxilla

Journal of Dentistry 2026



State of the Art Conference

See you at the Antwerp
Congress 2026!

Antwerp | October 16–17th, 2026



Digital
Dentistry
Society

ODIN WORKSHOP · STRASBOURG · 27 SEPTEMBER 2026

Thank you for your attention!

Current and Future Applications of XR and AI in Dentistry



Roberto Mangano

MSc Computer Engineering
Digital Dentistry Society Certified Speaker

www.digital-dentistry.org



Digital
Dentistry
Society